



Ophthalmic Stand SN-500

Serving Your Vision



Ophthalmic Stand SN-500

A versatile ophthalmic stand designed for reducing user's workload. A wide range of optical instruments can be set up according to user's preference.

Outstanding features

Overhead lamp

The overhead lamp knob on the control panel can dim up and down easily ensuring easy eye examination in the most suitable environment.

View Tester arm

A wide range of vertical movement and outstanding counterbalancing ensure smooth and easy operation with a view tester mounted.

Refractometer arm

A refractometer, etc. on an optional table top can move 220mm vertically.

Lowest arm

A slit lamp on an optional table top can move up and down easily by the up/down button under the table top.

* A power-supply voltage from each arm's power outlet is 50VA.

Control box

The center and left rechargeable instruments wells can accommodate a rechargeable ophthalmoscope each (DC 4V)*.

The right well can accommodate an ophthalmoscope connected to the binding post.

The output voltage of the binding post can be selected from 6 levels between AC 2.5V and AC 12V by the BD output voltage selection switch on the back.

* Some ophthalmoscopes cannot be charged depending on these battery types or sizes.

* Ophthalmoscopes with lithium ion batteries cannot be charged in the wells.

Specifications

Model name	SN-500 Ophthalmic Stand
Space required	535 × 535 mm
Weight	137 kg




Operation panel



SN-500 combined with Ophthalmic Chair UN-15

* VT-5 View tester, CP-40 Chart projector, UN-15 Ophthalmic chair, U08-04 Projector arm, and U12-05 Slit lamp table top are separately available.

 To ensure safe and correct handling please read the user manual before using.

•Description and appearance as detailed in this brochure may be subject to change as improvements are made to products.

•Colours as they appear in brochure photographs may be slightly different to actual products due to lighting conditions when photographs were taken, or print colours.



For The Americas, Asia-Pacific & Middle East

Takagi Seiko Co., Ltd

330-2 Iwafune, Nakano-shi, Nagano-ken, 383-8585, Japan
TEL : +81(0)269-22-6094 URL : <https://www.takagi-j.com>



For Europe & Africa

Takagi Ophthalmic Instruments Europe Ltd

Citylabs 1.0, Nelson Street, Manchester, M13 9NQ, UK
TEL : +44 (0)161 273 6330 URL : <https://www.takagieurope.com>



Your local distributor